

A KATENA BRAND

blink^omedical

ophthalmic single use instrument catalog



www.meditools.no

post@meditools.no

Tlf: 40 17 40 16

Blink Medical, et Katena merke

For mer informasjon om Blink Medical informasjon om alle Katena-produkter finner du online på www.katena.com. Ta kontakt med oss for å få full tilgang til alt undervisningsmaterieell og produktinformasjon. full tilgang til alt undervisningsmaterieell, produktinformasjone fra vårt direktesalgsstyrke og kundeserviceteam. For mer informasjon om Blink Medical, eller for å bestille: Informasjon om alle Katena-produkter finner du online på www.katena.com.

link Medical, a Katena brand

Blink Medical is a leading provider of sterile, single use ophthalmic instrumentation, instrument sets, and packs. Founded in 2004, Blink Medical has focused on providing high quality, cost effective single use instrumentation. Blink Medical was acquired by Katena Products, Inc. in 2018..

Blink Medical is a leading provider of sterile, single use ophthalmicInstrumentation, instrument sets, and packs. Founded in 2004, Blink Medical has focused on providing high quality, cost effective single use instrumentation. Blink Medical was acquired by Katena Products, Inc. in 2018.

For more information on Blink Medical, or to place an

Information on all Katena Products brands may be found online at www.katena.com.

Table of contents	
Product	Page
Speculums	4-5
Forceps	6-13
• Capsulorrhexis Forceps	6-7
• Tissue Forceps	8-9
• Tying Forceps	10
• Utility Forceps	11-12
• Jewelers Forceps	13
Manipulators & Choppers	14-17
Needle Holders	18
Calipers & Markers	19
Scissors	20-23
Muscle Hooks & Probes	24
Clamps	25
Curettes & Dilators	26
Other Instruments	27
Procedure Sets	28-29

Alle Blink Medical-produkter leveres i sterile i peel- pakker, trygt forpakket for sikker levering. Antall instrumenter per boks varierer avhengig kartpngent. Vennligst sjekk katalogen eller Katena.com for Æ fæÆå fræ fåp bekretet- det for å bekrefte antallet per boks.

Alle Blink Medical-produkter er beregnet for og anbefalles kun for engangsbruk. Sterilitetsholdbarhetn for e Blink Medical-produkter er tre (3) år fra steriliseringsdatoen. Se etiketten på esken for å finne utløpsdatoen for de forskjellige produktene.

hver boks med produkter.

Alle engangsinstrumentenes dimensjonale mål er å anse som circa mål. Sammenligbare

Sammenlignbare Katena produktnumre er vist er de nærmeste gjenbrukbare elementene og kan variere noe frenengangsinstrumentet.calecurely boxed for safe delivery.

Units per box varies depending upon product.

Speculums |

PEDIATRIC SPECULUM 30-3002

Spring with solid blade
Overall length 35 mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K1-5300
Sauer Infant Speculum



BARRAQUER SPECULUM 30-4000

Solid blade, 14 mm, angled, temporal
Overall length 40mm
Supplied sterile in boxes of 10
3 year shelf life


COMPARABLE KATENA REUSABLE: K1-5030
Barraquer Speculum Solid Large



KRATZ-BARRAQUER SPECULUM 31-1006

Wire blade, 15mm, angled, temporal
Overall length 40mm
Supplied sterile in boxes of 10
3 year shelf life


COMPARABLE KATENA REUSABLE: K1-5050
Kratz-Barraquer Speculum Large



BARRAQUER SPECULUM 31-2003

Wire blade, 15 mm, angled, temporal
Overall length 40 mm
Width 30 mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K1-5010
Barraquer Wire Speculum Large



LEIBINGER SPECULUM 31-3000

Solid blade, 15mm
Angled, temporal
Overall length 65mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K1-5678
Lieberman Speculum Solid Large



LEIBINGER SPECULUM 31-3002

Wire blade, 14 mm, angled, temporal
Overall length 65 mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K1-5675
Lieberman Speculum V-Wire



NEO-NATAL SPECULUM 31-3100

Wire blade, 6 mm
Angled, temporal
Overall length 34 mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K1-5350
Infant Wire Speculum 4mm Blade




Capsulorrhexis Forceps |

UTRATA CAPSULORRHEXIS FORCEPS

11-5012

Round handle, direct action
45° straight shaft
0.4mm sharp tip
Overall length 103 mm
Supplied sterile in boxes of 10
3 year shelf life




COMPARABLE KATENA REUSABLE: 08-01415
Utrata/Kershner Capsulorrhexis Forceps W/Cystotome Tips

UTRATA CAPSULORRHEXIS FORCEPS

11-5013

Round handle, direct action
Curved shaft
0.4mm sharp tip
Overall length 103 mm
Supplied sterile in boxes of 10
3 year shelf life




COMPARABLE KATENA REUSABLE: K5-5086
Masket Capsulorrhexis Forceps Round Handle

INAMURA CAPSULORRHEXIS FORCEPS

11-5020

Round handle, cross action
Curved shaft
0.4mm sharp tip
Overall length 105mm
Supplied sterile in boxes of 10
3 year shelf life




COMPARABLE KATENA REUSABLE: K5-5070
Inamura Capsulorrhexis Forceps Short

INAMURA CAPSULORRHEXIS FORCEPS

11-5021

Flat handle, cross action
Curved shaft
0.4mm sharp tip
Overall length 105mm
Supplied sterile in boxes of 10
3 year shelf




COMPARABLE KATENA REUSABLE: K5-5072
Inamura Capsulorrhexis Forceps Long

UTRATA CAPSULORRHEXIS FORCEPS

11-5022

Flat handle, direct action
45° straight shaft
0.4mm sharp tip
Overall length 105mm
Supplied sterile in boxes of 10
3 year shelf life




COMPARABLE KATENA REUSABLE: K5-5083
Nevyas Capsulorrhexis Forceps

MICS CAPSULORRHESIS FORCEPS

11-5023

- Micro-incision, 1.8mm
- Round handle
- Curved shaft, direct action
- 0.4mm sharp tip
- Overall length 97 mm
- Supplied sterile in boxes of 10
- 3 year shelf life



COMPARABLE KATENA REUSABLE: K5-5090
Giannetti MICS Capsulorrhexis Forceps Round Handle

UTRATA CAPSULORRHESIS FORCEPS

11-5024

- Flat handle
- Direct action, curved shaft
- 0.4mm sharp tip
- Overall length 105 mm
- Supplied sterile in boxes of 10
- 3 year shelf life



COMPARABLE KATENA REUSABLE: K5-5084
Masket Capsulorrhexis Forceps

FULCHER CAPSULORRHESIS FORCEPS

11-5025

- Flat handle
- Direct action, curved shaft
- 0.4mm rounded tip
- Overall length 105 mm
- Supplied sterile in boxes of 10
- 3 year shelf life



COMPARABLE KATENA REUSABLE: K5-5084
Masket Capsulorrhexis Forceps

UTRATA CAPSULORRHESIS FORCEPS

11-5026

- Flat handle
- Direct action, 45° shaft
- Overall length 105 mm
- Supplied sterile in boxes of 10
- 3 year shelf life





COMPARABLE KATENA REUSABLE: K5-5081
Utrata Capsulorrhexis Forceps

MICS CAPSULORRHESIS FORCEPS

11-5027

- Round handle
- 1.8mm incision, direct action
- 45° shaft
- Overall length 105mm
- Supplied sterile in boxes of 10
- 3 year shelf life





COMPARABLE KATENA REUSABLE: K5-5090
Giannetti MICS Capsulorrhexis Forceps Round Handle

Tissue Forceps |

BONN TISSUE FORCEPS

10-1105

1x2 teeth 0.25mm
6mm tying platform
Overall length 84mm
Supplied sterile in boxes of 10
3 year shelf life





COMPARABLE KATENA REUSABLE: K5-2000
Bonn Forceps Short 0.12mm

BONN MICRO TISSUE FORCEPS

10-1106

1x2 teeth 0.25mm
6mm tying platform
Overall length 106 mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K5-2020
Bonn Forceps Long 0.12mm

ST. MARTINS TISSUE FORCEPS

10-1203

1x2 teeth 0.6 mm
Overall length 90 mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K5-3562
Manhattan Eye & Ear Forceps

ST. MARTINS TISSUE FORCEPS

10-1204

1x2 teeth 0.3mm
Overall length 90mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K5-3562
Manhattan Eye & Ear Forceps

LISTER UTILITY FORCEPS

10-1302


1x2 teeth

6mm tying platform

Overall length 114mm

Supplied sterile in boxes of 10

3 year shelf life



COMPARABLE KATENA REUSABLE: K5-4300

Lester Fixation Forceps 1x2

HOSKINS MICRO FORCEPS

10-8009


0.3mm notch

5mm tying platform

Overall length 108mm

Supplied sterile in boxes of 10

3 year shelf life



COMPARABLE KATENA REUSABLE: K5-1718


Hoskin #18 Forceps Straight

COLIBRI TOOTHED FORCEPS

11-7010

1x2 teeth 0.25mm

Overall length 65 mm Supplied
sterile in boxes of 10 3 year shelf
life



COMPARABLE KATENA REUSABLE: K5-1520


Harms-Colibri Forceps

COLIBRI TOOTHED FORCEPS

11-7011

1x2 teeth 0.25mm

Overall length 110 mm Supplied
sterile in boxes of 10 3 year shelf
life



COMPARABLE KATENA REUSABLE: K5-1902


Harms-Colibri Forceps Long

Tying Forceps |

MCPHERSON TYING FORCEPS

10-9007

Straight plain tip
6mm tying platform
Overall length 85mm
Supplied sterile in boxes of 10
3 year shelf life




COMPARABLE KATENA REUSABLE: K5-5200
Jaffe Tying Forceps Straight

MCPHERSON TYING FORCEPS

10-9008

Straight plain tip
6mm tying platform
Overall length 110mm
Supplied sterile in boxes of 10
3 year shelf life




COMPARABLE KATENA REUSABLE: K5-5000
McPherson Tying Forceps Straight

MAX FINE TYING FORCEPS

11-6009

Curved, smooth tip
No platform
Overall length 110mm
Supplied sterile in boxes of 10
3 year shelf life




COMPARABLE KATENA REUSABLE: K5-5550
Fine Tying Forceps

KELMAN - MCPHERSON TYING FORCEPS

11-6010

8mm tying platform, angled
Overall length 85mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K5-5030
Kelman-McPherson Forceps Angled

Utility Forceps |

EPILATION FORCEPS

10-1300

Angled tips
 Overall length 91mm
 Supplied sterile in boxes of 10
 3 year shelf life



COMPARABLE KATENA REUSABLE: K5-6020
 Katena Cilia Forceps Slanted

MOORFIELDS SUTURE & UTILITY FORCEPS

10-1311

14mm serrated tying platform
 Overall length 107mm
 Supplied sterile in boxes of 10
 3 year shelf life



COMPARABLE KATENA REUSABLE: K5-4150
 Adson Dressing Forceps Serrated Delicate

MOSQUITO FORCEPS

10-1006

Ring handle
 Straight tip
 Overall length 100mm
 Supplied sterile in boxes of 10
 3 year shelf life




Picture is smaller than actual size

COMPARABLE KATENA REUSABLE: K5-9900
 Hartmann Mosquito Forceps Straight

MOSQUITO FORCEPS

11-1005

Ring handle
 Curved tip
 Overall length 100mm
 Supplied sterile in boxes of 10
 3 year shelf life




Picture is smaller than actual size

COMPARABLE KATENA REUSABLE: K5-9910
 Hartmann Mosquito Forceps Curved

Utility Forceps |

ADSON TOOTHED FORCEPS HR121

Dissecting 0.25mm 1X2 teeth
Overall length 120mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K5-4172
Adson Tissue Forceps 1x2 Teeth Standard

ADSON NON-TOOTHED FORCEPS HR125

Overall length 120mm
Serrated tips
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K5-4152
Adson Dressing Forceps Serrated Standard

ADSON-BROWN FORCEPS HR126

Overall length 120mm
Supplied sterile in boxes of 10
3 year shelf life



JEWELERS FORCEPS NO. 1

10-2001

Overall length 120mm
Standard, straight
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K5-6510
Jewelers Forceps #1 Straight



JEWELERS FORCEPS NO. 3

10-2003

Overall length 120mm
Delicate, straight
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K5-6530
Jewelers Forceps #3 Straight



JEWELERS FORCEPS NO. 5

10-2005

Overall length 115mm
Extra delicate, straight
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE:: K5-6550
Jewelers Forceps #5 Straight



JEWELERS FORCEPS NO. 5 ANGLED

11-2005

Overall length 112mm
45° angled shafts
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K5-6560



JEWELERS FORCEPS NO. 7

11-2007

Overall length 117.5mm
Standard, curved
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K5-6570
Jewelers Forceps #7 Curved



Manipulators/Choppers |

SINSKEY LENS HOOK

40-5002

Straight shaft

0.5mm tip

Overall length 124mm

Supplied sterile in boxes of 10

3 year shelf life



COMPARABLE KATENA REUSABLE: K3-5232

Sinskey II Lens Hook Straight Round

LESTER IOL MANIPULATOR

41-1006

Tip angled 45°

Angled at 12mm from tip

1mm branch

Overall length 125mm

Supplied sterile in boxes of 10

3 year shelf life



COMPARABLE KATENA REUSABLE: K3-2690

Lester IOL Manipulator Angled

Y-SHAPED NUCLEUS MANIPULATOR

41-2006

On vertical axis

Tip angled at 45°

Angled at 13mm from tip

Overall length 124mm

Supplied sterile in boxes of 10

3 year shelf life



COMPARABLE KATENA REUSABLE: K3-2400

Bechert Nucleus Rotator 7mm

NEUHANN SPATULA

41-3005

Vaulted spatula angled 11mm from tip

Supplied sterile in boxes of 10

3 year shelf life



COMPARABLE KATENA REUSABLE: K3-2388

Neuhann Nucleus Divider

AMBIDEXTROUS CHOPPER

41-4004

Tip angled at 45°

Angled at 11mm from tip

Overall length 115mm

Supplied sterile in boxes of 10

3 year shelf life



COMPARABLE KATENA REUSABLE: K3-2329

Nagahara Phaco Chopper

SADIQ CHOPPER

41-4005

Tip angled at 45°

Angled at 11mm from tip

Overall length 121mm

Supplied sterile in boxes of 10

3 year shelf life



COMPARABLE KATENA REUSABLE: K3-2396

Rosen Phaco Splitter 60 Degrees

CLAW CHOPPER

41-7009

Tip angled at 45°

Angled at 11mm from tip

Distal tip length 1mm

Overall length 120mm

Supplied sterile in boxes of 10

3 year shelf life



COMPARABLE KATENA REUSABLE: 08-14507

Steinert II Claw Chopper

Manipulators/Choppers |

MUSHROOM-TIPPED MANIPULATOR

41-8009

Tip angled at 45°

Angled at 11mm from tip

Overall length 123mm

Supplied sterile in boxes of 10

3 year shelf life



COMPARABLE KATENA REUSABLE: K3-2695
Hunkeler Lens Hook Hemi-Ball

SINSKEY/MUSHROOM-TIP DOUBLE ENDED

41-9010

Tips angled at 45°

Sinskey tip length 0.65mm

Overall length 123mm

Supplied sterile in boxes of 10

3 year shelf life



AYCLIFFE/SADIQ DOUBLE ENDED CHOPPER

41-9011

Tips angled at 45°

Aycliffe tip: 0.5mm

Sadiq tip: 1.5mm

Overall length 120mm

Supplied sterile in boxes of 10

3 year shelf life



MUSHROOM/SADIQ DOUBLE ENDED CHOPPER

41-9012

Tips angled at 45°
Overall length: 143mm
Sadiq tip: 1.5mm
Supplied sterile in boxes of 10
3 year shelf life



SINSKEY LENS HOOK

41-5001

Tip angled at 45°
Tip length: 0.5mm
Overall length 120mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K3-5233
Sinskey II Lens Hook Angled Round

DRYSDALE PADDLE

41-6007

Tip angled at 45°
Angled at 11mm from tip
Spatula 1mm
Overall length 120mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K3-2418
Drysdale Nucleus Manipulator


Needle Holders |

CASTROVIEJO NEEDLE HOLDER

90-7001

Straight with lock
12mm jaws
Overall length 135mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K6-1560
Castroviejo Needle Holder Straight W/Lock



CASTROVIEJO NEEDLE HOLDER

90-7002

Straight without lock
12mm jaws
Overall length 135mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K6-1510
Castroviejo Needle Holder Straight




CASTROVIEJO NEEDLE HOLDER

91-7003

Delicate curve without lock
12mm jaws
Overall length 135mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K6-1530
Castroviejo Needle Holder Curved



TROUTMAN NEEDLE HOLDER

91-7004

Curved without lock
6mm jaws
Overall length 105mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K6-3250
Troutman Needle Holder Curved



CALIPER & RULE

90-1014

Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K3-9000
Castroviejo Caliper Straight



DOUBLE ENDED MARKER

90-1015

3.5mm & 4.0mm ends
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K3-9250
Braunstein Caliper 3.5/4.0mm



THORNTON RING

90-1020

Globe fixation ring
Multiple blunt teeth
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K3-6161
Fine-Thornton Fixation Ring 13mm

DEGREE GAUGE

91-6000

11.7mm inside diameter
Marked in 5° increments
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K3-7904
Henderson Degree Gauge



ALIGNMENT MARKER

91-6001

3, 6, 9 o'clock placement
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K3-7908
Henderson Alignment Marker



TORIC IOL MARKER

91-6002

Two radial marks extending to 11mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K3-7912
Henderson Toric IOL Marker



Scissors |

WESTCOTT SCISSORS

20-4006

Straight tips

Overall length 120mm

Supplied sterile in boxes of 10

3 year shelf life

COMPARABLE KATENA REUSABLE: K4-3100

Westcott Tenotomy Scissors Right



VANNAS SCISSORS

20-6003

Straight tips

Overall length 95 mm

Supplied sterile in boxes of 10

3 year shelf life

COMPARABLE KATENA REUSABLE: K4-5000

Vannas Scissors Straight



IRIS SCISSORS

20-7001

Straight, rounded tips Overall length 100 mm Supplied sterile in boxes of 10 3 year shelf life



IRIS SCISSORS

20-7002

Straight, sharp tips

Overall length 100mm

Supplied sterile in boxes of 10

3 year shelf life

COMPARABLE KATENA REUSABLE: K4-7200

Eye Scissors Delicate Straight



WESTCOTT SCISSORS

21-4008

Curved, sharp tips
Overall length 120mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K4-4100
Westcott Stitch Scissors



WESTCOTT SCISSORS

21-4010

Curved, rounded tips
Overall length 120mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K4-3100
Westcott Tenotomy Scissors Right




WESTCOTT SCISSORS

21-4011

Dolphin nose tips
Overall length 120mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K4-3004
Westcott Tenotomy Scissors Medium



Scissors |

WESTCOTT SCISSORS SMALL

21-5001

Curved, sharp tips
Short blades
Overall length 120mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K4-4000
McPherson-Westcott Stitch Scissors

ONG SCISSORS

21-3010

Angled, sharp tips
Overall length 92mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K4-5020
Vannas Scissors Angled Forward

VANNAS SCISSORS

21-6004

Curved, sharp tips
Overall length 95mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K4-5010
Vannas Scissors Curved

IRIS SCISSORS

HR210

Curved, sharp tips
Overall length 100mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K4-7410
Eye Scissors Standard Curved

STRABISMUS SCISSORS

HR240

Curved, rounded tips
Overall length 115mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K4-9010
Knapp Strabismus Scissors Curved



STRABISMUS SCISSORS

HR241

Straight, rounded tips
Overall length 115mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K4-9000
Knapp Strabismus Scissors Straight



DRESSING STITCH SCISSORS

HR251

Sharp/Blunt tips
Overall length 115mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K4-9060
Utility Scissors



Probes & Muscle Hooks |

GRAEFE STRABISMUS HOOK 90-5004

Tip length 10mm
Overall length 140mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K3-6720
Graefe Muscle Hook #2

ALEXANDA MUSCLE HOOK 90-5005


Tip length 11mm
Overall length 130mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K3-6850
Green Muscle Hook

LIEBRICHT PROBE 90-2008


0000&000
Overall length 125mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K7-2510
Bowman Lacrimal Probe #0000-000

LIEBRICHT PROBE 90-2009

00&0
Overall length 125mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K7-2520
Bowman Lacrimal Probe #00-0

LIEBRICHT PROBE 90-2010

1&2
Overall length 125mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K7-2530
Bowman Lacrimal Probe #1-2

LIEBRICHT PROBE 90-2011

3&4
Overall length 125mm
Supplied sterile in boxes of 10
3 year shelf life



COMPARABLE KATENA REUSABLE: K7-2540
Bowman Lacrimal Probe #3-4

CHALAZION CLAMP

90-4005

Small

Overall length 90mm

Clamp diameter 11mm X 10mm

Supplied sterile in boxes of 10

3 year shelf life



COMPARABLE KATENA REUSABLE: K5-9500

Hunt Chalazion Forceps

CHALAZION CLAMP

90-4006

Medium

Overall length 90mm

Clamp diameter 18mm x 11mm

Supplied sterile in boxes of 10

3 year shelf life



COMPARABLE KATENA REUSABLE: K5-9600

Francis Chalazion Forceps

CHALAZION CLAMP

90-4007

Large

Overall length 90mm

Clamp diameter 23mm x 15mm

Supplied sterile in boxes of 10

3 year shelf life



COMPARABLE KATENA REUSABLE: K5-9650

Desmarres Chalazion Forceps

DESMARRES CLAMP

90-4008

Entropion left

Clamp diameter 24mm x 16mm

Overall length 95mm

Supplied sterile in boxes of 10

3 year shelf life



DESMARRES CLAMP

90-4009

Entropion right

Clamp diameter 24mm x 16mm

Overall length 95mm

Supplied sterile in boxes of 10

3 year shelf life



BULLDOG CLAMP

90-8012

Overall length 53mm

Supplied sterile in boxes of 10

3 year shelf life



COMPARABLE KATENA REUSABLE: K5-9850

Serrefine Small Straight

Curettes & Dilators |

CURETTE

90-3007

Double ended
2mm & 3mm cup
Overall length 130mm
Supplied sterile in boxes of 10
3 year shelf life



CURETTE

90-3008

Double ended
2mm & 5mm cup
Overall length 130mm
Supplied sterile in boxes of 10
3 year shelf life



CURETTE

90-3009

Double ended
2mm & 4mm cup
Overall length 130mm
Supplied sterile in boxes of 10
3 year shelf life



PUNCTAL DILATOR

90-1011

Double ended
Overall length 138mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K7-2300
Castroviejo Lacrimal Dilator



NETTLESHIPS DILATOR

90-1012

Overall length 125mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K7-2110
Wilder Lacrimal Dilator #1




EPITHELIAL SPATULA

51-0001

Spatula length 7mm
Overall length 113mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K2-3735
Fukasaku 'Hockey' Knife



ALCOHOL WELL

51-0002

Overall length 114mm
Inside diameter 9mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K3-1822
LASEK Alcohol Well 9.5mm



BARD-PARKER SCALPEL HANDLE #3

HR300

Overall length 135mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K20-1950
Bard-Parker Handle #3




BACKHAUS TOWEL CLIP

HR800

Overall length 95mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K5-9964
Backhaus Towel Clamp 3-1/2"



Picture is smaller than actual size

DESMARRES RETRACTOR

90-9013

14mm width
Overall length 130mm
Supplied sterile in boxes of 10
3 year shelf life

COMPARABLE KATENA REUSABLE: K1-8210
Desmarres Lid Retractor #1



CONJUNCTIVAL PROBE

CP01

For use with the SISTA™ technique.
(Small Incision Sub Tenon's Anesthesia)



SUB TENON'S SET 001A

Fine tip scissors

Contains:

1 x 21-4009 Westcott curved fine tip scissors


1 x 31-1006 Kratz-Barraquer angled temporal speculum

1 x Moorfields suture forceps

Protective wrap

Set supplied sterile in quantities of 8 per box

3 year shelf life



INTRAVITREORETINAL INJECTION SET 023

Contains:

1 x 30-1007 Kratz-Barraquer straight speculum


1 x 90-1015 Double ended scleral marker

1 x 10-1104 Bonn suture forceps

Protective wrap

Set supplied sterile in quantities of 8 per box

3 year shelf life



CHALAZION SET – SMALL CLAMP 092

1 x 90-4005 Chalazion clamp small

1 x 90-3009 Curette double ended med/small

No. 11 blade & handle

Protective wrap

Set supplied sterile in quantities of 8 per box

3 year shelf life



CHALAZION SET – MEDIUM CLAMP

111

Contains:

- 1 x 90-4006 Medium chalazion clamp
- 1 x 90-3009 Double ended medium/small curette
- No. 11 blade & handle
- Protective wrap
- Set supplied sterile in quantities of 8 per box
- 3 year shelf life



OPHTHALMIC SUTURE REMOVAL SET

210

Contains:

- 1 x 10-9007 McPherson tying forceps straight
- 1 x 21-6004 Vannas curved sharp tip scissors
- Protective wrap
- Set supplied sterile in quantities of 8 per box
- 3 year shelf life




CONJUNCTIVAL PROBE SUB TENON'S SET

270

Contains:

- 1 X 10-1301 Moorfields suture forceps
- 1 X 31-1006 Kratz-Barraquer angled temporal speculum
- 1 X CP01 conjunctival probe
- Protective wrap
- Set supplied sterile in quantities of 10 per box
- 3 Year shelf life



A KATENA BRAND

blink medical

MEDI-tOOLS AS
Vestvollveien 6j
2020 Skedsmokorset



www.katena.com