



VISCOT[®]

MEDICAL, LLC.

The *Leader* In Surgical Marker Innovation

Surgical Markers • Medication Labels • Neutral
Zones Drapes • Urinals

www.viscot.com

P:1.800.221.0658

F:973.887.3961

viscotcs@viscot.com

2022

Viscot is proud to be a veteran owned small business (VOSB)



post@meditools.no
Tlf: 40 17 40 16

v02-11.17.21

ABOUT US

Viscot Medical, LLC, is the global leader in surgical site marking innovations. Viscot is committed to providing the healthcare industry with high quality disposable products at an economical price. We are proud to serve our global customers as a U.S. manufacturer for over 40 years.

Viscot Medical is dedicated to listening to clinicians and coming up with site marking solutions to meet each clinical need.

Viscot Medical is a veteran-owned, small business.

ISO 13485 CERTIFIED

TABLE OF CONTENTS

click on page numbers for quick navigation

About Our Inks	3
About Our Tips	4
About Our Barrels	5
Traditional Ink Surgical Markers	7
XL Prep Resistant Ink Markers	15
Non-Surgical Aesthetic Skin Markers	19
Specialized Site Markers	21
Utility Markers	25
Pre-Printed Medication Labels	28
Sharps & Staff Safety	32
Additional Products	34
Catalog Marker Lookup	37



TRADITIONAL INK

Traditional ink is Gentian Violet ink and has been used in the industry to mark patients for over 40 years. It can withstand betadine and iodine scrubs. AllSkin[®] (pg 10) features a brilliant blue option available for visibility on all skin tones.

Ideal for procedures above the neck or where betadine and iodine are used.

Featured on pages 8-13



GV FREE INK

Viscot Medical's GV Free ink is a gentian violet free alternative to Viscot's Traditional ink. Eliminates concerns over Gentian Violet while maintaining the classic purple color.

Ideal for procedures above the neck or where betadine and iodine are used.

Featured on page 14




XL PREP RESISTANT INK

Viscot XL Prep Resistant Ink is visible up to 10x longer vs. traditional gentian violet ink. Able to withstand a 60 second prep of ChloroPrep, XL Ink remains visible against all common alcohol based prep solutions. Available in gentian violet free red and green ink for enhanced visibility on all skin tones.

Ideal for any procedure where aggressive preps are used.

Featured on pages 16-18



EZ REMOVABLE INK™

Viscot Medical's EZ Removeable Ink is designed to come off easily with an alcohol swab or even warm soap and water. Available in red, green, and white.

Ideal for non-surgical aesthetic procedures.

Featured on page 20



WATERPROOF UTILITY INK

Black Utility ink is designed for writing on medication labels. Smudge-proof, waterproof, and fluid resistant.

Ideal for writing on labels, glass, vinyl and other surfaces.

Featured on page 26-27



ULTRA FINE TIPS

Viscot Medical's soft rounded Ultra Fine tips create fine lines while eliminating concerns over causing trauma to delicate tissues and reducing frustration from sharp marker points snagging making it difficult to draw.

Ideal for facial procedures such as blepharoplasty procedures.



REGULAR TIPS

This is the most common tip style found on the market. Designed to provide versatility and a middle-ground between fine lines and thick lines.

Ideal for most general procedures and for correct site marking.



BOLD TIPS

Viscot Medical's bold tips cover a wide surface area and create broad lines. Additional tipsurface makes it more difficult for the marker to clog and appear dry.

Ideal for abdominoplasty procedures and vascular procedures.



BROAD TIPS

Viscot Medical's broad tip is unique to our 1404 and 1405 markers. This tip was specifically designed to replace gentian violet ink stamp pads. The paintbrush like tip coats instruments easily and completely.

Ideal for coating surgical instruments with ink but especially cataract procedures.



MINI MARKERS

Viscot Medical's Mini Markers are economical and eco-friendly. They're about half the size of a standard marker so they're about half the cost and produce half the waste.

Perfect for promoting single patient use.



PRECISION BARREL

Our standard barrels designed with a comfortable grip in mind. Precision Barrels feel and perform like a fountain pen. Featured on our Precision Markers and full size XL markers.



WIDE BARREL

Thicker barrels that feel more like a standard marker and less like a pen. All Wide Barrel Markers have a clip cap.



VALUE BARREL

This barrel feels like a pencil when gripped. Value barrels are thinner and shorter than Precision and Wide barrels, but not as small as a Mini Marker. Less plastic means more cost savings.



MULTI-TIP BARRELS

Features two tips one on each end. Multi-Tip barrels differ depending on the exact product numbers but are an economical, convenient and eco-friendly alternative by providing two markers in one.

Usually featuring different tips, the unique Blephmarker[®] (pg 9) features two ultra fine tips, offering a fresh tip per eyelid.



PRESSURIZED BARREL *

Found in Viscot Medical's Visflow Markers. Designed to keep marking even through fluids on wet tissue.

Ideal for abdominoplasty procedures, breast procedures, and liposuction.

*patent pending



TRADITIONAL INK SURGICAL MARKERS

Pre-Op & Intraoperative Markers Filled
with Gentian Violet Ink or GVFree alternatives.
Ideal for procedures where Betadine or
Iodine are used.



NEW

TOR-MARK-10-5

Disposable Cataract Marking Kit

Mini Max Regular Tip Traditional Ink Marker, Dial with Level and Indexed Cap 5/cs





1400/1401 PRECISION SURGICAL MARKER
Regular Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1400-100	Sterile Regular Tip w/ Ruler	100
1400-50L	Sterile Regular Tip w/ 4 Labels & Ruler	50
1400-100L	Sterile Regular Tip w/ 4 Labels & Ruler	100
1401-100	Non-Sterile Regular Tip	100



1402/1403 PRECISION SURGICAL MARKER
Ultra Fine Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1402-50	Sterile Ultra Fine Tip w/ Ruler	50
1402-100	Sterile Ultra Fine Tip w/ Ruler	100
1403-100	Non-Sterile Ultra Fine Tip	100



1404/1405 PRECISION SURGICAL MARKER
Broad Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1404-50	Sterile Broad Tip w/ Ruler	50
1404-100	Sterile Broad Tip w/ Ruler	100
1405-100	Non-Sterile Broad Tip	100



1406/1407 PRECISION SURGICAL MARKER / ALLSKIN
Brilliant Blue Ink | Regular Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1406-100	Non-Sterile Regular Tip	100
1407-50	Sterile Regular Tip w/ Ruler	50
1407-100	Sterile Regular Tip w/ Ruler	100



1422 DUAL TIP BARREL
Regular Tip + Traditional Ink | Fine Tip + Permanent Utility Ink

PART #	DESCRIPTION	CASE QTY
1422-100	Non-Sterile 1 Regular & 1 Fine Tip	100
1422S-100	Sterile 1 Regular & 1 Fine Tip	100
1422SR-50	Sterile w/ Ruler 1 Regular & 1 Fine Tip	50
1422SR-100	Sterile w/ Ruler 1 Regular & 1 Fine Tip	100
1422SRL9-50	Sterile 1 Regular & 1 Fine Tip w/ Ruler & 9 Blank Labels	50
1422SRL9-100	Sterile 1 Regular & 1 Fine Tip w/ Ruler & 9 Blank Labels	100



1424 TWIN TIP BARREL
Blephmarker: 2 Ultra Fine Tips | Traditional Gentian Violet Ink

PART #	DESCRIPTION	CASE QTY
1424-100	Non-Sterile 2 Ultra Fine Tips	100
1424SR-100	Sterile 2 Ultra Fine Tips w/ Ruler	100



1425 TWIN TIP BARREL

Regular Tip & Ultra Fine Tip | Traditional Gentian Violet Ink

PART #	DESCRIPTION	CASE QTY
1425-100	Non-Sterile 1 Regular & 1 Ultra Fine Tip Sterile 1	100
1425SR-100	Regular & 1 Ultra Fine Tip w/ Ruler	100



1436 VALUE SURGICAL MARKER

Ultra Fine Tip | Value Barrel

PART #	DESCRIPTION	CASE QTY
1436-100	Non-Sterile Ultra Fine Tip	100
1436SR-50	Sterile Ultra Fine Tip w/ Ruler	50
1436SR-100	Sterile Ultra Fine Tip w/ Ruler	100
1436SRL9-100	Sterile Ultra Fine Tip w/ Ruler & 9 Blank Labels	100



1437 VALUE SURGICAL MARKER

Regular Tip | Value Barrel

PART #	DESCRIPTION	CASE QTY
1437-100	Non-Sterile Regular Tip	100
1437SR-50	Sterile Regular Tip w/ Ruler	50
1437SR-100	Sterile Regular Tip w/ Ruler	100
1437SRL9-50	Sterile Regular Tip w/ Ruler & 9 Blank Labels	50
1437SRL9-100	Sterile Regular Tip w/ Ruler & 9 Blank Labels	100



1438 VALUE SURGICAL MARKER
Regular Tip | Value Barrel | Tail on Cap

PART #	DESCRIPTION	CASE QTY
1438SR-100	Sterile Regular Tip, Clip Cap w/ Ruler	100



1451 MINI SURGICAL MARKER
Regular Tip | Mini Barrel

PART #	DESCRIPTION	CASE QTY
1451-200	Non-Sterile Regular Tip - 2 Plastic Containers of 100	200
1451-1000	Non-Sterile Regular Tip - 10 Plastic Containers of 100	1,000
1451SR-100	Sterile Regular Tip w/ Ruler	100
1451W-200	Non-Sterile Regular Tip Individually Wrapped (Proof of Use) - 4 Plastic Containers of 50 Markers	200
1451W-1000	Non-Sterile Regular Tip Individually Wrapped (Proof of Use) - 20 Plastic Containers of 50 Markers	1,000



1460 STANDARD TRADITIONAL INK SURGICAL MARKER
Regular Tip | Wide Barrel

PART #	DESCRIPTION	CASE QTY
1460-100	Non-Sterile Regular Tip	100
1460SRL-100	Sterile Regular Tip w/ Ruler & Labels	100
1460SRL9-100	Sterile Regular Tip w/ Ruler & 9 Labels	100
		12.



1480 STANDARD SURGICAL MARKER
Regular Tip | Wide Barrel

PART #	DESCRIPTION	CASE QTY
1480-100	Non-Sterile Regular Tip	100



VF400 VISFLOW PRESSURIZED SKIN MARKER *
Bold Tip | Pressurized Barrel

PART #	DESCRIPTION	CASE QTY
VF400SR-25	Sterile Bold Tip w/ Ruler	25
VF400SR-50	Sterile Bold Tip w/ Ruler	50



VF500 VISFLOW ULTRA SQUEEZABLE INK MARKER *
Regular Tip | Pressurized Barrel

PART #	DESCRIPTION	CASE QTY
VF500SR-50	Sterile Regular Tip w/ Ruler	50
VF500SR-100	Sterile Regular Tip w/ Ruler - 2/Bx of 50 Markers	100

*patent pending

A Gentian Violet Free Alternative To Traditional
Gentian Violet Free Markers. Complies with California Prop 65.



1487 VALUE GVFREE INK SKIN MARKER
Regular Tip | Value Barrel

PART #	DESCRIPTION	CASE QTY
1487GVFSR-100	Sterile Regular Tip w/ Ruler	100
1487GVFSRL-100	Sterile Regular Tip w/ Ruler & 9 Blank Labels	100



1481 MINI GVFREE INK SKIN MARKER
Regular Tip | Mini Barrel

PART #	DESCRIPTION	CASE QTY
1481GVF-200	Non-Sterile Regular Tip - 2 Plastic Containers of 100	200
1481GVF-1000	Non-Sterile Regular Tip - 10 Plastic Containers of 100	1,000
1481GVFSR-100	Sterile Regular Tip w/ Ruler	100

VLSCOT[®]
MEDICAL, LLC.

XL PREP RESISTANT INK MARKERS

Designed For Pre-Operative Site Marking
Where Aggressive Preps Are Used.
Visible Throughout Prepping with All Prep
Including Chlorhexidine & Alcohol Based Preps.





1447 XL PREP RESISTANT INK SKIN MARKER

Regular Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1447XL-50	Non-Sterile Regular Tip - 1 Plastic Container of 50	50
1447XL-50	Non-Sterile Regular Tip - 1 Plastic Container: 18 Red, 16 Green, & 16 Purple	50
1447XLSR-100	Sterile Regular Tip w/ Ruler	100



1448 XL PREP RESISTANT INK SKIN MARKER

Ultra Fine Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1448XL-50	Non-Sterile Ultra Fine Tip - 1 Plastic Container of 50	50
1448XLSR-100	Sterile Ultra Fine Tip w/ Ruler	100



1449 XL PRE-SURGICAL MARKER

Bold Tip | Precision Barrel

PART #	DESCRIPTION	CASE QTY
1449XL-50	Non-Sterile Bold Tip - 1 Plastic Container of 50	50
1449XLSR-100	Sterile Bold Tip w/ Ruler	100



1450 MINI MARKER XL PREP RESISTANT INK SKIN MARKER
Regular Tip | Mini Barrel

PART #	DESCRIPTION	CASE QTY
1450XL-200	Non-Sterile Regular Tip - 2 Plastic Containers of 100	200
1450XL-1000	Non-Sterile Regular Tip - 10 Plastic Containers of 100	1,000
1450XLSR-100	Sterile Regular Tip w/ Ruler	100
1450XLW-200	Non-Sterile Regular Tip Individually Wrapped (Proof of Use) - 4 Plastic Containers of 50	200
1450XLW-1000	Non-Sterile Regular Tip Individually Wrapped (Proof of Use) - 20 Plastic Containers of 50	1,000



1456 MINI MARKER XL PREP RESISTANT INK SKIN MARKER
Ultra Fine Tip | Mini Barrel

PART #	DESCRIPTION	CASE QTY
1456XL-200	Non-Sterile Ultra Fine Tip - 2 Plastic Containers of 100	200
1456XL-1000	Non-Sterile Ultra Fine Tip - 10 Plastic Containers of 100	1,000
1456XLSR-100	Sterile Ultra Fine Tip w/ Ruler	100
1456XLW-200	Non-Sterile Ultra Fine Tip Individually Wrapped (Proof of Use) - 4 Plastic Containers of 50	200



1447 XL PREP RESISTANT RED INK SKIN MARKER
Regular Tip | Precision Barrel | Gentian Violet Free

PART #	DESCRIPTION	CASE QTY
1447XLR-50	Non-Sterile Regular Tip - 1 Plastic Container	50
1447XLMC-50	Non-Sterile Regular Tip - 1 Plastic Container: 18 Red, 16 Green, & 16 Purple	50



1447 XL PREP RESISTANT GREEN INK SKIN MARKER

Regular Tip | Precision Barrel | Gentian Violet Free

PART #	DESCRIPTION	CASE QTY
1447XLG-50	Non-Sterile Regular Tip - 1 Plastic Container	50
1447XLMC-50	Non-Sterile Regular Tip - 1 Plastic Container: 18 Red, 16 Green, & 16 Purple	50
1447XLGSR-100	Sterile Regular Tip w/ Ruler	100



1470 XL PREP RESISTANT RED INK MINI SKIN MARKER

Regular Tip | Mini Barrel | Gentian Violet Free

PART #	DESCRIPTION	CASE QTY
1470XLRED-1000	Non-Sterile Regular Tip - 10 Plastic Containers of 100	1,000
1470XLREDSR-100	Sterile Regular Tip w/ Ruler	100
1470XLRED-200	Non-Sterile Regular Tip - 2 Plastic Containers of 100	200



1477 XL PREP RESISTANT RED INK VALUE SKIN MARKER

Regular Tip | Value Barrel | Gentian Violet Free

PART #	DESCRIPTION	CASE QTY
1477XLREDSR-100	Sterile Regular Tip w/ Ruler	100
1477XLREDSRL9-100	Sterile Regular Tip w/ Ruler & 9 Blank Labels	100

VLSCOT[®]
MEDICAL, LLC.

NON-SURGICAL AESTHETIC SKIN MARKERS

Marking Non-Surgical Aesthetic Procedures
Like Consultations, Botox, & Fillers.





1440 EZ REMOVABLE WHITE INK SKIN MINI MARKER

Regular Tip | Mini Barrel | White Ink

PART #	DESCRIPTION	CASE QTY
1440-30	Non-Sterile Regular Tip - 1 Plastic Container	30
1440-240	Non-Sterile Regular Tip - 8 Plastic Containers of 30	240



1444 EZ REMOVABLE GREEN INK SKIN MINI MARKER

Regular Tip | Mini Barrel | Green Ink

PART #	DESCRIPTION	CASE QTY
1444-30	Non-Sterile Regular Tip - 1 Plastic Container	30
1444-240	Non-Sterile Regular Tip - 8 Plastic Containers of 30	240



1445 EZ REMOVABLE RED INK SKIN MINI MARKER

Regular Tip | Mini Barrel | Red Ink

PART #	DESCRIPTION	CASE QTY
1445-30	Non-Sterile Regular Tip - 1 Plastic Container	30
1445-240	Non-Sterile Regular Tip - 8 Plastic Containers of 30	240

VLSCOT[®]
MEDICAL, LLC.

SPECIALIZED SITE MARKERS

Alternatives To Surgical Markers.





AM106 ACCU-MARK RADIOPAQUE SITE MARKER
Marking Under Fluoroscopy

PART #	DESCRIPTION	CASE QTY
AM106-10	Non-Sterile w/ Traditional Surgical Ink	10
AM106S-10	Sterile w/ Traditional Surgical Ink	10
AM106-50	Non-Sterile w/ Traditional Surgical Ink	50
AM106S-50	Sterile w/ Traditional Surgical Ink	50
AM106XL-50	Non-Sterile w/ XL Prep Resistant Ink	50
AM106-100	Non-Sterile w/ Traditional Surgical Ink	100
AM106-250	Non-Sterile w/ Traditional Surgical Ink	250



X SITE "YES" RADIOPAQUE STICKERS
Endourologic Procedural Marking


PART #	DESCRIPTION	CASE QTY
XS100	Non-Sterile	100
XS250	Non-Sterile	250



LASER TAT MARKERS

For Ophthalmic Laser Site Marking

PART #	DESCRIPTION	CASE QTY
TEL-400	Non-Sterile "Laser" in Red w/ White Background 1" Diameter - 100 Per Box / 4 Boxes / 1 Cs	400



EYE IMPLANT TAT MARKERS

For Cataract Implant Site Marking

PART #	DESCRIPTION	CASE QTY
TEY-400	Non-Sterile Ophthalmic - 100 Per Box / 4 Boxes / 1 Cs	400



RADIOLOGY SITE MARKERS

For Oncology and Interventional Radiology Site Marking

PART #	DESCRIPTION	CASE QTY
TXR-400	Non-Sterile Radiation Oncology - 100 Per Box / 4 Boxes / 1 Cs	400



RIGHT AND LEFT SITE MARKERS

For Ophthalmic, Pediatric and Oncology Site Marking

PART #	DESCRIPTION	CASE QTY
TLR5050	Non-Sterile 50 Green "Left" & 50 Blue "Right" -100 Per Box / 4 Boxes / 1 Cs	400



SMALL 1/2" YES SITE MARKERS

Cataract, Pediatric & Oncology Site Marking

PART #	DESCRIPTION	CASE QTY
TYS-100W	Non-Sterile Temporary 1/2" in Green "Yes" - 100 Markers w/ 100 Pre-Moistened Wipes	100
TYS-400	Non-Sterile Temporary 1/2" in Green "Yes" - 100 Per Box / 4 Boxes / 1 Cs	400



LARGE 1" YES SITE MARKERS



Cataract, Pediatric & Oncology Site Marking

PART #	DESCRIPTION	CASE QTY
TYL-100W	Non-Sterile Temporary 1" in Green "Yes" - 100 Markers w/ 100 Pre-Moistened Wipes	100
TYL-400	Non-Sterile Temporary 1" in Green "Yes" - 100 Per Box / 4 Boxes / 1 Cs	400



UTILITY MARKERS

Marking Medication Labels & Surfaces
Such as Plastic, Glass, Metal, etc.





1408 WATERPROOF PERMANENT MARKER
Regular Tip | Wide Barrel

PART #	DESCRIPTION	CASE QTY
1408-50	Non-Sterile Regular Tip	50
1408SRL-100	Sterile Regular Tip w/ Ruler & 4 Blank Labels	100



1411 WATERPROOF PERMANENT MARKER
Fine Tip | Wide Barrel

PART #	DESCRIPTION	CASE QTY
1411-100	Non-Sterile Fine Tip	100
1411S-100	Sterile Fine Tip	100
1411SRL9-50	Sterile Fine Tip w/ Ruler & 9 Blank Labels	50
1411SRL9-100	Sterile Fine Tip w/ Ruler & 9 Blank Labels	100



1416 WATERPROOF PERMANENT MARKER
Ultra Fine Tip | Wide Barrel

PART #	DESCRIPTION	CASE QTY
1416-100	Non-Sterile Ultra Fine Tip	100
1416S-100	Sterile Ultra Fine Tip	100



1422 DUAL TIP SURGICAL INK MARKER & PERMANENT INK MARKER

Regular Tip + Traditional Ink | Fine Tip + Permanent Utility Ink

PART #	DESCRIPTION	CASE QTY
1422-100	Non-Sterile Regular Tip & Non-Sterile Fine Tip	100
1422S-100	Sterile Regular Tip & Fine Tip	100
1422SR-50	Sterile Regular Tip & Fine Tip w/ Ruler	50
1422SR-100	Sterile Regular Tip & Fine Tip w/ Ruler	100
1422SRL9-50	Sterile Regular Tip & Fine Tip w/ Ruler & 9 Blank Labels	50
1422SRL9-100	Sterile Regular Tip & Fine Tip w/ Ruler & 9 Blank Labels	100



1453 WATERPROOF PERMANENT MARKING MINI MARKER



Regular Tip | Mini Barrel

PART #	DESCRIPTION	CASE QTY
1453-100	Non-Sterile Regular Tip	100
1453SRL9-100	Sterile Regular Tip w/ Ruler & 9 Blank Labels	100



PRE-PRINTED MEDICATION LABELS

Custom labels available for no additional cost.



Ancef	Lidocaine _____%
Ancef	Lidocaine _____%
Adenosine _____mcg/ml	Nipride _____mcg/ml
Adenosine _____mcg/ml	Nitro _____mcg/ml
Contrast Name _____mgI/ml	Normal Saline 0.9%
Contrast Name _____mgI/ml	Normal Saline 0.9%
Heparin _____units/ml	Verapamil _____mcg/ml
Heparin _____units/ml	Verapamil _____mcg/ml
Heparinized Saline _____units/ml	
Heparinized Saline _____units/ml	

CATH LAB COLOR MEDICATION LABEL KIT

PART #	DESCRIPTION	CASE QTY
VCATH-20	Sterile 20 Color Kit	100
VCATH-22R-20	Sterile 20 Color Kit w/ 1422 Dual Tip Marker & Ruler	100

AIR	Mitomycin _____mg/ mL
Antibiotic	Moxifloxacin 0.5% (Vigamox)
BSS	Perfluoron
BSS	SF ₆ _____%
C ₃ F ₈ _____%	Sterile Water
Decadron _____mg/mL	Triamcinolone PF _____mg/mL
Dexamethasone _____mg/mL	Trypan Blue
Lidocaine _____% MPF	Vitrise _____units/mL
Miochol	
Miostat	

OPHTHALMIC COLOR MEDICATION LABEL KIT

PART #	DESCRIPTION	CASE QTY
VEYE-20	Sterile 20 Color Kit	100
VEYE-22R-20	Sterile 20 Color Kit w/ 1422 Dual Tip Marker & Ruler	100

Contrast Name _____ mgI/ml	Bupivacaine _____ %
Contrast Name _____ mgI/ml	Bupivacaine _____ %
Heparin _____ units/ml	Lidocaine 1% w/ Sodium Bicarb Amt _____ Exp _____
Heparin _____ units/ml	Lidocaine 1% w/ Sodium Bicarb Amt _____ Exp _____
Heparinized Saline _____ units/ml	Normal Saline 0.9%
Heparinized Saline _____ units/ml	Normal Saline 0.9%
Lidocaine _____ %	1% Lidocaine HCL w/ Epinephrine 1:100,000 Amt _____ Exp _____
Lidocaine _____ %	1% Lidocaine HCL w/ Epinephrine 1:100,000 Amt _____ Exp _____
Nipride _____ mcg/ml	
Nitro _____ mcg/ml	

**INTERVENTIONAL RADIOLOGY COLOR
MEDICATION LABEL KIT**

PART #	DESCRIPTION	CASE QTY
VIR-20	Sterile 20 Color Kit	100
VIR-22R-20	Sterile 20 Color Kit w/ 1422 Dual Tip Marker & Ruler	100

Bupivacaine _____ %	Depomedrol _____ mg/ml
Bupivacaine _____ %	Depomedrol _____ mg/ml
Lidocaine _____ %	Betamethasone _____ mg/ml
Lidocaine _____ %	Betamethasone _____ mg/ml
Contrast Name _____ mgI/ml	Dexamethasone _____ mg/ml
Contrast Name _____ mgI/ml	Dexamethasone _____ mg/ml
Kenalog _____ mg/ml	Verapamil _____ mcg/ml
Normal Saline 0.9%	Verapamil _____ mcg/ml
Normal Saline 0.9%	
Triamcinolone _____ mg/mL	

**PAIN MANAGEMENT COLOR
MEDICATION LABEL KIT**

PART #	DESCRIPTION	CASE QTY
VPAIN-20	Sterile 20 Color Kit	100
VPAIN-22R-20	Sterile 20 Color Kit w/ 1422 Dual Tip Marker & Ruler	100

Acetic Acid	Methergine
Bupivacaine _____% with Epinephrine	Methylene Blue 50mg/5ml
Bupivacaine _____%	Mineral Oil
Estradiol 0.01% Cream	Monssel's Solution
Hydrogen Peroxide 3%	Normal Saline 0.9%
Iodine Scrub	Silver Sulfadizine 1% (Silvadene)
Iodine Prep	Vasopressin 20u/ml
KY Jelly	WATER
Lidocaine 1% Plain	
Lugol's Solution	

OBGYN COLOR MEDICATION LABEL KIT

PART #	DESCRIPTION	CASE QTY
VOBGYN-20	Sterile 20 Color Kit	100
VOBGYN-22R-20	Sterile 20 Color Kit w/ 1422 Dual Tip Marker & Ruler	100



BLANK MEDICATION LABELS

PART #	DESCRIPTION	CASE QTY
V20063DB-DTR	Sterile 10 Blank Labels w/ 1422 Dual Tip Marker & Ruler	100

VLSCOT[®]
MEDICAL, LLC.

SHARPS & STAFF SAFETY

Reusable - Autoclavable





3" x 9" SOFFZONE NEUTRAL ZONE MAT *

PART #

FNZ309RED

DESCRIPTION

3"x9" Red Tray w/ Raised Contours for Sharps

CASE QTY

6



6" x 9" SOFFZONE NEUTRAL ZONE MAT *

PART #


FNZ609RED

DESCRIPTION

6"x9" Red Mat w/ Raised Contours for Sharps

CASE QTY

6



SOFFZONE SYRINGE TRAY *

PART #

FSY309RED

DESCRIPTION

3"x9" Red Tray w/ Tabs for Holding Syringes

CASE QTY

6

*patent pending

VLSCOT[®]
MEDICAL, LLC.

ADDITIONAL PRODUCTS





MILLIE FEMALE URINAL

PART #	DESCRIPTION	CASE QTY
201000-12	Female Urinal	12
201036	Female Urinal	36



**UNIVERSAL URINAL
w/ Lid**

PART #	DESCRIPTION	CASE QTY
201000-L-12	Universal Urinal w/ Leak Resistant Lid - Individually Bagged w/ Instructions	12
201036-L	Universal Urinal w/ Leak Resistant Lid - Individually Bagged w/ Instructions	36



MALE URINAL

PART #	DESCRIPTION	CASE QTY
221000	Male Urinal	50
		35.

Rulers & labels included in sterile packaging:

Viscot Medical Sterile markers come individually packaged often with a fluid resistant ruler and sometimes with back up labels. Rulers are flexible fluid resistant coated paper and include both metric and U.S. measurements units (15 cm / 6 in).

Blank back up labels included with markers come in quantities of 4 or 9. They are paper based labels and are designed as a back up to our standard pre-printed medication labels.



1409 METRIC ONLY (15CM) RULER

PART #	DESCRIPTION	CASE QTY
1409-1000	Non-Sterile	1,000

1410 METRIC/INCH (15CM/6IN) RULER

PART #	DESCRIPTION	CASE QTY
1410-100	Non-Sterile	100
1410-1000	Non-Sterile	1,000
1410S-300	Sterile	300

R24 24IN FOLDED RULER

PART #	DESCRIPTION	CASE QTY
R24-1000	Non-Sterile	1,000

PRODUCT #	1400/1401	1402/1403	1404/1405	1406/1407	1408	1411	1416	1422	1424	1425	1436	1437	1438	1440	1444	1445	1447XL	1447XLG	1447XLMC	1447XLR	1448XL	1449XL	1450XL	1451	1453	1456XL	1460	1470XLRED	1477XLRED	1480	1481GVF	1487GVF	TOR-MARK-10-5	VF400	VF500					
PAGE #	9	9	9	10	26	26	26	10/27	10	11	11	11	12	20	20	20	16	18	16/17/18	17	16	16	17	12	27	17	12	18	18	13	14	14	8	13	13					
TRADITIONAL INK (GENTIAN VIOLET)	█	█	█					█	█	█	█	█	█											█										█	█	█				
GVFREE (GENTIAN VIOLET FREE TRADITIONAL INK)																																								
ALLSKIN (BRILLIANT BLUE)				█																																				
XL PREP RESISTANT INK																	█	█	█	█	█	█	█			█														
XL PREP RESISTANT INK - COLORS (GENTIAN VIOLET FREE)																		█	█	█	█	█	█					█												
EZ REMOVABLE INK															█	█																								
UTILITY INK (LABEL MARKING)					█	█	█	█																		█														
MINI														█	█												█													
FULL SIZE PRECISION	█	█	█			█												█	█	█	█	█	█																	
FULL SIZE WIDE					█		█																				█													
FULL SIZE VALUE											█	█	█																█											
FULL SIZE DOUBLE TIP								█	█																															
PRESSURIZED																																						█	█	
ULTRA FINE		█						█	█	█	█	█	█	█	█	█					█						█													
REGULAR	█			█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█	█
BOLD																																							█	█
BROAD			█																																					

www.viscot.com